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The emergence of a new series specialized for UHPLC!!

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
■ Provide total solution in 1 series

The Develosil UHPLC series is a series specialized only for UHPLC adopting 1.6 μ m particles. For the lineup, we have five types of C30, C18, C8, C1, HILIC, you can get all the selectivity with this one series and make it possible to achieve total separation.

For ordinary analysis centered on hydrophobic compounds, C18 is the first choice and scouting to C8 and C1 according to the required retention. And when targeting a highly polar compound or a highly ionic compound, selecting C30 and C1 will exhibit different characteristics from general purpose C18. Also, select HILIC to thoroughly hold high polar compounds eluting at the front! It improves the efficiency of separation analysis.

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
■ The selected silica gel base material


The Develosil UHPLC series is not a conventional silica gel base material, but a highly pure silica gel base material is selected. This selection further reduces unnecessary metal coordination and others to produce a sharp peak shape, and the reinforced silica gel base material shows high durability even under high pressure conditions peculiar to UHPLC analysis, ensuring stable data at all times To do.

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
■ The selected silica gel base material

The five types of Develosil UHPLC columns control the bond density to take advantage of individual characteristics with new and new functions. In particular, since it is possible to analyze by increasing the organic solvent, it is possible to improve the sensitivity at the time of using LC/MS and to withdraw from the use of ion pair reagent in peptide analysis, etc., resulting in revolutionary results You can get.

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

The retention obtained by the degree of hydrophobicity can easily control the analysis cycle by changing the length of the alkyl chain.

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Changing the pH can easily change the pattern of retention and separation. In addition, establishment of analytical method that does not use ion pair reagent greatly improves column load reduction and cleaning efficiency.

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
■Bold method development

In the Develosil UHPLC series, in addition to the inversion mode, the "Develosil UHPLC HILIC" of the HILIC mode was lineup. HILIC is widely used as a method for efficiently analyzing highly polar compounds by raising the composition of organic solvents, but by selecting 1.6 μ m particles, analysis in a shorter time is possible.

● Develosil UHPLC HILICによる塩基性医薬品成分の分析 Analysis of basic pharmaceutical ingredients by Develosil UHPLC HILIC

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Further information can be obtained from the website.

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Develosil UHPLC 1.6 μ m Metal-free Series

Size	C30	C18	C8	C1	HILIC
2.0x35mm	201-I20035MFW	202-I20035MFW	203-I20035MFW	204-I20035MFW	205-I20035MFW
2.0x50mm	201-I20050MFW	202-I20050MFW	203-I20050MFW	204-I20050MFW	205-I20050MFW
2.0x75mm	201-I20075MFW	202-I20075MFW	203-I20075MFW	204-I20075MFW	205-I20075MFW
2.0x100mm	201-I20100MFW	202-I20100MFW	203-I20100MFW	204-I20100MFW	205-I20100MFW
2.0x150mm	201-I20150MFW	202-I20150MFW	203-I20150MFW	204-I20150MFW	205-I20150MFW

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〒489-0004 愛知県瀬戸市日の出町15
 TEL:0561-48-1853 FAX:0561-48-1434
 E-mail: info@develosil.net
 Web: www.develosil.net

NOMURA CHEMICAL CO., LTD.

15, Hinode-cho, Seto, 489-0004, Japan
 TEL: +81-561-48-1853 FAX: +81-561-48-1434
 E-mail: info@develosil.net
 Web: www.develosil.net/en/

Develosil USA

Develosil USA

10060 Carroll Canyon Rd. Ste. 100 San Diego, CA 92131
 Phone: 858-800-2433
 Web: https://develosil.us/