

Technical Report No.32

Analysis of peptides/proteins and monoclonal antibodies by UHPLC part II

1. Column for polymer compounds

Proteins and monoclonal antibodies are compounds with molecular weights of tens of thousands to hundreds of thousands. A column with a pore size of about 300 Å is required to analyze these.

The FlexFire WP series has a lineup of 4 types at 300 Å and 1 type as a column dedicated to monoclonal antibodies (1000 Å), for a total of 5 types, and you can select according to your purpose.

In this report, we will introduce the usefulness of wide pore columns, targeting higher molecular compounds than the previous one.

2. Specification


	FlexFire mAb-RP	FlexFire WP C4	FlexFire WP C18	FlexFire WP C8	FlexFire WP C1
Particle size	2.6µm, 5µm	2.6µm, 5µm	2.6µm, 5µm	2.6µm, 5µm	2.6µm, 5µm
Chemistry	Butyl	Butyl	Octadecyl	Octyl	Trimethyl
Surface area	24m ² /g	170m ² /g	170m ² /g	170m ² /g	170m ² /g
Pore volume	1.4mL/g	1.4mL/g	1.4mL/g	1.4mL/g	1.4mL/g
Pore diameter	100nm	30nm	30nm	30nm	30nm
Carbon	1.3%	5%	15%	7%	3%
End-cap	O	O	O	O	O
pH	pH1-10	pH1-10	pH1-10	pH1-10	pH1-10
Temperature	~80°C	~80°C	~80°C	~80°C	~80°C
Pressure range	2.6µm: 600bar (=60Mpa=8,702psi)	2.6µm: 600bar (=60Mpa=8,702psi)	2.6µm: 600bar (=60Mpa=8,702psi)	2.6µm: 600bar (=60Mpa=8,702psi)	2.6µm: 600bar (=60Mpa=8,702psi)
	5µm: 300bar (=30Mpa=4,351psi)	5µm: 300bar (=30Mpa=4,351psi)	5µm: 300bar (=30Mpa=4,351psi)	5µm: 300bar (=30Mpa=4,351psi)	5µm: 300bar (=30Mpa=4,351psi)

3. mAb recovery rate


The recovery rate of mAb (150KDa) in each column was calculated.

It shows a high recovery rate in every column.

Column	R ²	Recovery (%)
FlexFire WP C1	0.99886	104
FlexFire WP C4	0.99961	103
FlexFire WP C8	0.99795	89.8
FlexFire WP C18	0.99823	98.1
FlexFire mAb-RP	0.99886	101



4. Comparison of separation between myoglobin and impurities.



Conditions;

Column: FlexFire WP Series, 2.6 μ m (2.0x50mm)
 Mobile phase: A) Water+0.1%TFA B) Acetonitrile+0.1%TFA

Gradient:

min	mL/min	%A	%B	Curve
0.00	0.3	80	20	
8.40	0.3	40	60	6
8.42	0.3	80	20	6

Temperature: 40 $^{\circ}$ C

Detection: UV280nm

Sample: Myoglobin

Injection volume: 0.2 μ L


System: Waters ACQUITY UPLC H-Class PLUS

Mixer: 100 μ L

The separation of myoglobin and impurities is good with either column.

Since high resolution can be obtained with FlexFire mAb-RP, WP C8, and WP C4, it can be expected to be effective in separating multiple components.

5. Analysis of Lactalbumin



Conditions;

Column: FlexFire C18, 2.6 μ m (2.0x50mm)
 FlexFire WP C4, 2.6 μ m (2.0x50mm)
 FlexFire mAb-RP, 2.6 μ m (2.0x50mm)

Mobile phase: A) Water + 0.1%TFA B) Acetonitrile + 0.1%TFA

Gradient:

min	mL/min	%A	%B	Curve
0.00	0.3	80	20	
8.40	0.3	40	60	6
8.42	0.3	80	20	6

Temperature: 40 $^{\circ}$ C


Detection: UV280nm

Sample: 1. α -Lactalbumin (0.34mg/mL)
 2. β -Lactoglobulin (1.00mg/mL)

Injection volume: 2.0 μ L

System: Waters ACQUITY UPLC H-Class PLUS

Mixer: 100 μ L




Compared to columns with a pore size of 120 \AA , columns with 300 \AA and 1000 \AA have better peak separation and sensitivity.

Especially at 1000 \AA , clearer peak separation is achieved.



6. Separation and comparison of multi-component proteins.



Conditions:

Column: FlexFire mAb-RP, 2.6 μm (2.0x50mm)
 Mobile phase: A) Water + 0.1%TFA B) Acetonitrile + 0.1%TFA

Gradient:	min	mL/min	%A	%B	Curve
	0.00	0.3	80	20	
	8.40	0.3	40	60	6
	8.42	0.3	80	20	6

Temperature: 40°C
 Detection: UV210nm
 Sample:
 1. Ribonuclease A (13.7KDa)
 2. Cytochrome C (12.4KDa)
 3. Lysozyme (14.3KDa)
 4. BSA (66.3KDa)
 5. Myoglobin (11.2KDa)
 6. Catalase (220KDa)

Injection volume: 2.0 μL

System: Waters ACQUITY UPLC H-Class PLUS
 Mixer: 100 μL



Separation with 10KDa-60KDa is possible even with a column with a pore diameter of 300 Å.

It can be seen that among the 200KDa protein, the column with the larger pore size achieved the separation predominantly.

By using FlexFire mAb-RP, not only separation between macromolecules but also detection of impurities contained in proteins can be expected.



7. Analysis of mAb



Conditions;

Column: FlexFire mAb-RP, 2.6 μ m (2.0x50mm)
 FlexFire WP C4, 2.6 μ m (2.0x50mm)
 Mobile phase: A) Water + 0.1%TFA B) Acetonitrile + 0.1%TFA


Gradient:	min	mL/min	%A	%B	Curve
	0.00	0.3	80	20	
	8.40	0.3	40	60	6
	8.42	0.3	80	20	6

Temperature: 40 $^{\circ}$ C
 Detection: UV280nm
 Sample: NISTmAb
 Injection volume: 1.0 μ L

System: Waters ACQUITY UPLC H-Class PLUS
 Mixer: 100 μ L

Technical Report No. 30 focused on applications with a pore diameter of 300 \AA , but this report uses a column with a pore diameter of 1000 \AA . A 1000 \AA column can find clearer peaks and gain new insights.

Technical Report No. 30 has an example of analysis using a 300 \AA column.



Conditions;

Column: FlexFire mAb-RP, 2.6 μ m (2.0x50mm)
 Mobile phase: A) Water + 0.1%TFA B) Acetonitrile + 0.1%TFA

Gradient:	min	mL/min	%A	%B	Curve
	0.00	0.3	80	20	
	8.40	0.3	40	60	6
	8.42	0.3	80	20	6

Temperature: 40 $^{\circ}$ C
 Detection: UV280nm
 Sample: NISTmAb, Reduced NISTmAb
 Injection volume: 1.0 μ L

System: Waters ACQUITY UPLC H-Class PLUS
 Mixer: 100 μ L


Analysis of mAb reducer

20 mM DTT was mixed 1: 1 with NISTmAb and reacted at 40 degree-2 hr.

Technical Report No. 30 has an example of analysis using a 300 \AA column.



8. Peptide mapping (BSA digest)



Conditions;

Column: FlexFire C18, 2.6µm (2.0x50mm)
 FlexFire C18, 1.6µm (2.0x50mm)
 FlexFire WP C18, 2.6µm (2.0x50mm)
 FlexFire WP C18, 1.6µm (2.0x50mm)
 FlexFire WP C4, 2.6µm (2.0x50mm)
 FlexFire mAb-RP, 2.6µm (2.0x50mm)

Mobile phase: A) Water + 0.1%TFA B) Acetonitrile + 0.1%TFA

Gradient:

min	mL/min	%A	%B	Curve
0.00	0.3	80	20	
8.40	0.3	40	60	6
8.42	0.3	80	20	6

Temperature: 40°C

Detection: UV280nm

Sample: BSA Digest

Injection volume: 1.0µL

System: Waters ACQUITY UPLC H-Class PLUS


Mixer: 100µL

The Thermo Scientific™ SMART Digest™ Kit was used to digest the protein.

When including intact analysis, it is necessary to select a column suitable for the molecular weight, but if only peptides that have been reduced in molecular weight by digestion, a standard column such as 120 Å can be selected.

※ Thermo Scientific™ SMART Digest™ Kit is a registered product of Thermo Fisher Scientific Co., Ltd. .

8-1. mAb digest



Conditions;

Column: FlexFire mAb-RP, 2.6µm (2.0x50mm)

Mobile phase: A) Water + 0.1%TFA B) Acetonitrile + 0.1%TFA

Gradient:

min	mL/min	%A	%B	Curve
0.00	0.3	80	20	
8.40	0.3	40	60	6
8.42	0.3	80	20	6

Temperature: 80°C

Detection: UV280nm

Sample: NISTmAb Digest

Injection volume: 1.0µL

System: Waters ACQUITY UPLC H-Class PLUS

Mixer: 100µL


mAb is a compound with a molecular weight of approximately 150 KDa.

The usefulness of intact mAb 1000Å has been found in the analysis of "7. Monoclonal antibody".

Even in peptide mapping, you can complete the process from intact analysis to digestion analysis with a single column without changing the column.



8-2. Hemoglobin digest



Conditions;
 Column: FlexFire mAb-RP, 2.6µm (2.0x50mm)
 Mobile phase: A) Water + 0.1%TFA B) Acetonitrile + 0.1%TFA


min	mL/min	%A	%B	Curve
0.00	0.3	80	20	
8.40	0.3	40	60	6
8.42	0.3	80	20	6

Gradient:

Temperature: 40°C
 Detection: UV280nm
 Sample: Hemoglobin Digest
 Injection volume: 1.0µL

System: Waters ACQUITY UPLC H-Class PLUS
 Mixer: 100µL

8-3. ADH digest



Conditions;
 Column: FlexFire mAb-RP, 2.6µm (2.0x50mm)
 Mobile phase: A) Water + 0.1%TFA B) Acetonitrile + 0.1%TFA


min	mL/min	%A	%B	Curve
0.00	0.3	80	20	
8.40	0.3	40	60	6
8.42	0.3	80	20	6

Gradient:

Temperature: 40°C
 Detection: UV280nm
 Sample: ADH Digest
 Injection volume: 1.0µL

System: Waters ACQUITY UPLC H-Class PLUS
 Mixer: 100µL

8-4. Insulin digest



Conditions;
 Column: FlexFire mAb-RP, 2.6µm (2.0x50mm)
 Mobile phase: A) Water + 0.1%TFA B) Acetonitrile + 0.1%TFA

min	mL/min	%A	%B	Curve
0.00	0.3	80	20	
8.40	0.3	40	60	6
8.42	0.3	80	20	6


Gradient:

Temperature: 40°C
 Detection: UV280nm
 Sample: Insulin Digest
 Injection volume: 1.0µL

System: Waters ACQUITY UPLC H-Class PLUS
 Mixer: 100µL



9. Analysis of synthetic RNA





Conditions:
 Column: FlexFire C18, 1.6µm (2.0x50mm)
 FlexFire WP C18, 1.6µm (2.0x50mm)
 Mobile phase: A) 100mM HFIP + 10mM TEA B) Methanol
 Gradient:

min	mL/min	%A	%B	Curve
0.00	0.3	95	5	
8.60	0.3	80	20	6
8.61	0.3	95	5	6

Temperature: 50°C
 Sample: RNA
 Injection volume: 0.2µL
 System: Waters ACQUITY UPLC H-Class PLUS
 Mixer: 100µL

By using a wide pore column, it can be detected more sensitively than a 120Å column.

10. Analysis of pBR322 digest



Conditions:
 Column: FlexFire C18, 1.6µm (2.0x50mm)
 FlexFire WP C18, 1.6µm (2.0x50mm)
 Mobile phase: A) 0.1M TEAA, pH7.0 B) 0.1M TEAA, pH7.0/ACN=80/20
 Gradient:

min	mL/min	%A	%B	Curve
0.00	0.2	42.5	57.5	
20.0	0.2	15.5	84.5	6
20.1	0.2	42.5	57.5	6

Temperature: 50°C
 Sample: pBr322 MspI Digest
 Injection volume: 10µL
 System: Waters ACQUITY UPLC H-Class PLUS
 Mixer: 100µL

In the analysis of polymer compounds, analysis that cannot be achieved with FlexFire C18 (120 Å) can be performed by using 300 Å "FlexFire WP C18".

■ お問い合わせ/Contact us



野村化学株式会社
〒489-0004 愛知県瀬戸市日の出町15
Tel: 0561-48-1853 Fax: 0561-48-1434
e-mail: info@develosil.net

Nomura Chemical Co., Ltd.
15, Hinode-cho, Seto, 489-0004, Japan
Tel: +81-561-48-1853 Fax: +81-561-48-1434
e-mail: info@develosil.net



Develosil USA
10060 Carroll Canyon Rd. Ste. 100 San Diego, CA 92131
Phone: 858-800-2433
Web: <https://develosil.us/>

