



**Concept**

## What's FlexFire ?



*FlexFire's Flex expresses Flexible, and Fire expresses the strength of fire. The logo was created with the motif of the phoenix as a symbol of the fire.*

*Forty years have passed since Nomura Chemical started "Develosil". We made a big decision in this memorable year. It was a renewal of silica gel. Making new silica gel is not easy. However, by reviewing all the know-how and user feedback that we have accumulated over the past 40 years, we were able to understand what was needed.*

*I hope that the FlexFire series born from this will set a new history for Nomura Chemical.*

## Get started with FlexFire! !

## *Transfer from the Develosil series*

The FlexFire series can be easily transferred from those who use the Develosil ODS series, HG / UG series, XG series, and HSR series to those who develop new methods. In addition, differences in the degree of separation that occur during transfer are supported by the method development technology that we have cultivated so far.

**ODS Series**

**HG, UG Series**

**XG Series**

**HSR**  
Develosil®

**UHPLC**  
Develosil®




## New Silica-gel


*What is the difference between the previous products and the FlexFire series?*

*-It differs in the silica gel selected as shown above.*

*Fig. 1 shows the results of comparing the performance of conventional silica gel and new silica gel. As shown in this figure, we were able to obtain silica gel with uniform particles and few impurities.*



**Fig.2 TEM image of silica gel for FlexFire (1.6 $\mu$ m-5 $\mu$ m)**



### Analytical conditions;

Column: FlexFire HILIC, 3 $\mu$ m (4.6x150mm)  
Develosil Silica 100-3, 3 $\mu$ m (4.6x150mm)  
Mobile phase: Acetonitrile/Water=90/10  
Flow rate: 1.0mL/min  
Temperature: 40°C  
Detection: CAD  
Sample: Allantoin  
Injection volume: 1.0 $\mu$ L

System: Thermo Fisher SCIENTIFIC UltiMate 3000

**Fig.1 Comparison of materials**

## System matching

*“Is your system HPLC? or UHPLC?” With the FlexFire series, the grain size can be selected according to the system. For any system, the FlexFire series is the most effective for your system.*

*Φ4.6mm and Φ1.0mm columns have been added to the FlexFire series to suit many systems.*

### *【 FlexFire series and system matching example 】*

| <i>Category</i> | <i>Compatible i.d.</i>           | <i>Compatible particle size</i> | <i>Verification system</i>  |
|-----------------|----------------------------------|---------------------------------|---|
| <i>HPLC</i>     | <i>Φ 4.6mm</i>                   | <i>2.6μm, 5μm</i>               | <i>Waters alliance</i>  |
| <i>UHPLC</i>    | <i>Φ 2.0mm</i><br><i>Φ 1.0mm</i> | <i>1.6μm, 2.6μm, 5μm</i>        | <i>Waters H-Class PLUS</i><br><i>Shimadzu Nexra X3</i><br><i>Agilent 1290 Infinity II</i><br><i>Thermo Vanquish H</i> |

*This table shows the results based on our own verification in normal use. Φ4.6mm can be used in UHPLC systems, but the flow cell and piping need to be replaced.*


## FlexFire's Line up

| Mode                          | Product name    | Particle size |       |     | Pore diameter<br>(nm) | Surface area<br>(m <sup>2</sup> /g) | Pore volume<br>(mL/g) | Carbon<br>(%) |
|-------------------------------|-----------------|---------------|-------|-----|-----------------------|-------------------------------------|-----------------------|---------------|
|                               |                 | 1.6µm         | 2.6µm | 5µm |                       |                                     |                       |               |
| Reversed phase                | FlexFire C18    | ○             | ○     | ○   | 11                    | 340                                 | 1.0                   | 22            |
|                               | FlexFire AQ C18 | ○             | ○     | ○   | 11                    | 340                                 | 1.0                   | 8.5           |
|                               | FlexFire C8     | ○             | ○     | ○   | 11                    | 340                                 | 1.0                   | 12            |
|                               | FlexFire C1     | ○             | ○     | ○   | 11                    | 340                                 | 1.0                   | 5.5           |
|                               | FlexFire C8     | ○             | ○     | ○   | 11                    | 340                                 | 1.0                   | 12            |
|                               | FlexFire C30    | ○             | ○     | ○   | 11                    | 340                                 | 1.0                   | 11            |
| Reversed phase<br>(Wide Pore) | FlexFire WP C18 |               | ○     | ○   | 30                    | 170                                 | 1.4                   | 15            |
|                               | FlexFire WP C8  |               | ○     | ○   | 30                    | 170                                 | 1.4                   | 7             |
|                               | FlexFire WP C4  |               | ○     | ○   | 30                    | 170                                 | 1.4                   | 5             |
|                               | FlexFire WP C1  |               | ○     | ○   | 30                    | 170                                 | 1.4                   | 3             |
|                               | FlexFire mAb-RP |               | ○     | ○   | 100                   | 24                                  | 0.8                   | 1.3           |
| SEC                           | FlexFire 120SEC |               | ○     | ○   | 11                    | 340                                 | 1.0                   | 6             |
|                               | FlexFire 300SEC |               | ○     | ○   | 30                    | 170                                 | 1.4                   | 9             |
| HILIC                         | FlexFire HILIC  | ○             | ○     | ○   | 11                    | 340                                 | 1.0                   | —             |

*This list is as of November 2021. The pH range has now been updated to pH 1-10 (except for HILIC). The FlexFire series will be further updated, such as WPC4, which targets high molecular compounds such as proteins and monoclonal antibodies.*

*And one of the big decisions is that the FlexFire series has only a 2.0mm ID column. This is a future-ready reform with the spread of UHPLC systems.*


## Durability and bond density



### Conditions;

Column: FlexFire C18, 1.6µm 2.0x50mm  
Size: Acetonitrile/10mM NH<sub>4</sub>OH, pH10.5=60/40  
Mobile phase: 0.5mL/min  
Flow rate: 40°C  
Temperature: UV254nm  
Detection:  
Sample: 1.Uracil (0.01mg/mL)  
2.Naphthalene (0.1mg/mL)  
Injection volume 0.16µL

**Fig.3 Retention time maintenance rate with respect to liquid flow time**




**Fig.4 Bond density of each chemistry**

*Silica-gel undergoes cleavage of the bonded phase under acidic conditions. And, the silica gel dissolves under alkaline conditions. However, the new silica-gel substrate of the FlexFire series has improved its strength. Figure 3 demonstrates that FlexFire is sufficiently durable at pH 10.5.*

*Until now, the durability changed according to the bond density, but tough silica gel has sufficient durability even for low-density chemistry.*


# UHPLC Method Transfer

The FlexFire series is available in particle sizes of 1.6 $\mu$ m, 2.6 $\mu$ m and 5 $\mu$ m. The particle size can be selected according to the system. And method transfer from 5 $\mu$ m to 1.6 $\mu$ m can be performed easily.



**Fig.5 Pressure and HETP vs Flow Rate for FlexFire C18**

**Fig.5 Image model of method transfer**




## Conditions;

Column: FlexFire C18 (2.0x50mm)  
Mobile phase: Acetonitrile/Water=60/40  
Flow rate: 0.2mL/min~1.0mL/min  
Temperature: 40 $^{\circ}$ C  
Detection: UV254nm  
Sample: Naphthalene  
Injection volume: 0.16 $\mu$ L

**Fig.6 Pressure and HETP vs Flow Rate for FlexFire C18**



## Example of method transfer



### Conditions:

Column FlexFire C18, 5µm (2.0x150mm)

FlexFire C18, 1.6µm (2.0x50mm)

Mobile phase: Acetonitrile/Water=60/40

Flow rate: 5µm: 0.189mL/min

1.6µm: 0.591mL/min

Temperature: 40°C

Detection: UV254nm

Sample: 1.Uracil

2.Naphthalene

Injection volume: 5µm: 0.3µL

1.6µm: 0.1µL

System: Waters ACQUITY UPLC H-Class PLUS


*The advantage of UHPLC method transfer is a significant reduction in analysis time.*

*In this example, about 8 minutes were reduced to 0.8 minutes by method transfer.*

*Moreover, method transfer can be easily performed with the attached software.*

**Fig.5 Example simple method transfer**

# Comparison of insulin peak shapes



**Fig.6 Comparison of Insulin peak shapes**

*Insulin was analyzed using a 0.1% formic acid mobile phase. The peak shape of insulin differs depending on the column. Insulin uses a TFA mobile phase to significantly improve the peak shape, but the introduction to LC / MS is severe.*

## Conditions:

Column: FlexFire C18, 1.6µm (2.0x50mm)

Mobile phase: A) Water + 0.1%HCOOH

B) Acetonitrile + 0.1%HCOOH

|           |      |        |    |    |       |
|-----------|------|--------|----|----|-------|
| Gradient: | min  | mL/min | %A | %B | Curve |
|           | 0.00 | 0.5    | 80 | 20 | 5     |
|           | 5.04 | 0.5    | 40 | 60 | 5     |
|           | 5.05 | 0.5    | 80 | 20 | 5     |

Temperature: 40°C

Detection: UV280nm


Sample: Insulin, Human, recombinant(0.97mg/mL)

Injection volume: 0.3µL

System: Thermo Fisher SCIENTIFIC Vanquish\_H

Mixer: 10µL


## Comparison of insulin peak shapes with Various mobile phases



**Fig.7 Comparison of Insulin peak shapes with various mobile phase**

*If good results are not obtained with 0.1% formic acid, it is desirable to switch to DFA or TFA. If ion suppression or contamination of the ion source affects, it is necessary to reduce the effect of the ion pair reagent by mixing with formic acid or the like.*

## Conditions that can be used with 100% aqueous mobile phase



**Fig.8 Stop-flow test under 100% aqueous mobile phase conditions**

### Conditions:


|                   |   |
|-------------------|---|
| Column:           | FlexFire AQ C18, 2.6μm  |
| Size:             | 2.0x100mm, Stainless  |
| Mobile Phase:     | 10mM HCOONH <sub>4</sub>  |
| Flow rate:        | 0.3mL/min   |
| Temperature:      | 40°C  |
| Detection:        | UV260nm   |
| Sample:           | 1.Cytosine (53μg/mL)<br>2.Uracil (50μg/mL)<br>3.Guanine (52μg/mL)<br>4.Thymine (50μg/mL)<br>5.Adenine (50μg/mL) |
| Injection volume: | 0.2μL   |

*The most important thing to consider with 100% aqueous mobile phase is reproducibility.*

*The mechanism is very simple. By stopping the pump, the mobile phase taken into the pores escapes. Then, even when the liquid sending is started again, the mobile phase is not taken into the pores.*

*This can be avoided by controlling the bond density. FlexFire AQ C18, C30, C1 controls the bond density so that it can be used in 100% aqueous system.*

## Overwhelming holding power with 100% aqueous mobile phase



*The 100% aqueous mobile phase has a significant effect on improving the separation of the front part.*

*In particular, the longer the alkyl chain, the stronger the retention of the hydrophobic compound.*


*That is, a column that can be used with an aqueous 100% mobile phase can separate low polar compounds from high polar compounds in a well-balanced manner.*

### Conditions;

|                   |   |
|-------------------|---|
| Column:           | FlexFire C30, 1.6µm (2.0x50mm)                  |
|                   | Other column, 2µm (2.0x50mm)                    |
|                   | ※Columns for analysis of highly polar compounds |
| Mobile phase:     | 25mM Ammonium phosphate, pH7.0                  |
| Flow rate:        | 0.5mL/min                                       |
| Temperature:      | 40℃   |
| Detection:        | UV260nm   |
| Sample:           | 1.dCTP 2.dTTP 3.dGTP 4.dATP                     |
| Injection volume: | 0.2µL   |

**Fig.9 Analysis of dNTP with 100% aqueous mobile phase**

## *S/N up with short alkyl chain*




**Fig.10 S/N comparison (C18vs C1)**

*Perhaps some people do not know the existence of C1. C1 is, so to speak, the last point of the reversed phase. Isn't C1 an image that is weak and unsuitable for separation? C1 does not have a short retention, and can be very good if it is considered to have low adsorption.*

*Now, does the data right here look like that? The retention and separation of the front part is not inferior to C18. An immediate understanding is the retention of hydrophobic compounds (Naphthalene). However, S / N has overwhelmingly high C1. And the peak that was not C18 can be confirmed in C1.*

*This has a very good effect on metabolic analysis. The combination of UHPLC methods and low adsorption further accelerates analysis time savings.*

## C8 column for sharp peaks



### Conditions;

Column: FlexFire C8, 2.6µm (2.0x50mm)  
 Develosil C8-UG, 3µm (2.0x50mm)  
 Other C8, 3µm (2.0x50mm)

Mobile phase: A) Water + 0.1%HCOOH  
 Flow rate: B) Acetonitrile + 0.1%HCOOH

### Gradient:

| min  | mL/min | %A | %B | Curve |
|------|--------|----|----|-------|
| 0.00 | 0.3    | 90 | 10 |       |
| 2.93 | 0.3    | 75 | 25 | 6     |
| 2.96 | 0.3    | 90 | 10 | 6     |

### Detection:

UV260nm

### Sample:

1. Gly-Tyr 2. Val-Tyr-Val 3. Angiotensin II  
 4. Met-Enk 5. Leu-Enk

### Injection volume:

0.2uL

### System:

Waters ACQUITY UPLC H-Class PLUS

**Fig.11 Separation comparison of Peptide**

*FlexFire C8 was developed for retention and separation. It is particularly excellent in separation among C8 columns of the same class. When the peak stands sharply, the S / N tends to be high. If the separation at C18 is insufficient, there is a good chance that it can be improved by changing to C8.*


## Metal-Free column

*Metal free column was added from FlexFire series. Compounds containing phosphorus groups adsorb to metals. Until now, the use of a phosphate buffer improved the peak shape, but it had a major problem.*


*It cannot be introduced into LC / MS.  
So we looked at metal-free columns and could get very good results with volatile buffers.*

*FlexFire's metal-free columns have a PEEK tube built into a stainless steel tube. This ensures high durability even under high pressure of UHPLC columns such as 1.6µm particle size.*

**Stainless**




**Metal-Free**






## Example using metal-free column (Flavin mononucleotide)



### Conditions;


Column: FlexFire AQ C18, 2.6 $\mu$ m (2.0x100mm) : Stainless  
FlexFire AQ C18, 2.6 $\mu$ m (2.0x100mm) : Metal Free  
Mobile phase: Acetonitrile/25mM HCOONH<sub>4</sub>=10/90  
Flow rate: 0.3mL/min  
Temperature: 40°C  
Detection: UV254nm  
Sample: Flavin mononucleotide (0.52mg/mL)  
Injection volume: 0.2 $\mu$ L

System: Waters ACQUITY UPLC H-Class PLUS



*Flavin mononucleotide*

## Example using metal-free column (AMP, ADP, ATP)



### Conditions:

Column: FlexFire AQ C18, 2.6 $\mu$ m (2.0x100mm) : Stainless  
FlexFire AQ C18, 2.6 $\mu$ m (2.0x100mm) : Metal Free

Mobile phase: 10mM HCOONH<sub>4</sub>

Flow rate: 0.3mL/min


Temperature: 40°C

Detection: UV260nm


Sample: 1. ATP (0.16mg/mL)  
2. ADP (0.17mg/mL)  
3. AMP (0.16mg/mL)

Injection volume: 0.2 $\mu$ L


System: Waters ACQUITY UPLC H-Class PLUS



AMP




ADP



ATP

## Example using metal-free column (dNTP)



Conditions:

Column: FlexFire AQ C18, 2.6 $\mu$ m (2.0x100mm) : Stainless  
 FlexFire AQ C18, 2.6 $\mu$ m (2.0x100mm) : Metal Free

Mobile phase: 10mM HCOONH<sub>4</sub>

Flow rate: 0.3mL/min

Temperature: 40°C

Detection: UV260nm

Sample: 1. dCTP (13 $\mu$ M)





2. dTTP (13 $\mu$ M)

3. dGTP (13 $\mu$ M)

4. dATP (13 $\mu$ M)

Injection volume: 0.2 $\mu$ L

System: Waters ACQUITY UPLC H-Class PLUS



# Analysis of biopolymer compounds

*In recent years, biopharmacy has grown into a very large market. In particular, nucleic acids and antibodies have grown remarkably, and with the spread of corona-virus, they have become familiar keywords to the general public.*


*In addition, HPLC / UHPLC is essential for these analyses, and manufacturers strive daily for better results.*

*We have succeeded in developing and commercializing a new wide pore column in the FlexFire series.*

## FlexFire series lineup

|                | FlexFire mAb-RP                    | FlexFire WP C4                     | FlexFire WP C18                    | FlexFire WP C8                     | FlexFire WP C1                     |
|----------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| Particle size  | 2.6µm, 5µm                         | 2.6µm, 5µm                         | 2.6µm, 5µm                         | 2.6µm, 5µm                         | 2.6µm, 5µm                         |
| Chemistry      | Butyl                              | Butyl                              | Octadecyl                          | Octyl                              | Trimethyl                          |
| Surface area   | 24m <sup>2</sup> /g                | 170m <sup>2</sup> /g               | 170m <sup>2</sup> /g               | 170m <sup>2</sup> /g               | 170m <sup>2</sup> /g               |
| Pore volume    | 1.4mL/g                            | 1.4mL/g                            | 1.4mL/g                            | 1.4mL/g                            | 1.4mL/g                            |
| Pore diameter  | 100nm                              | 30nm                               | 30nm                               | 30nm                               | 30nm                               |
| Carbon         | 1.3%                               | 5%                                 | 15%                                | 7%                                 | 3%                                 |
| End-cap        | O                                  | O                                  | O                                  | O                                  | O                                  |
| pH             | pH1-10                             | pH1-10                             | pH1-10                             | pH1-10                             | pH1-10                             |
| Temperature    | ~80°C                              | ~80°C                              | ~80°C                              | ~80°C                              | ~80°C                              |
| Pressure range | 2.6µm: 600bar<br>(=60Mpa=8,702psi) | 2.6µm: 600bar<br>(=60Mpa=8,702psi) | 2.6µm: 600bar<br>(=60Mpa=8,702psi) | 2.6µm: 600bar<br>(=60Mpa=8,702psi) | 2.6µm: 600bar<br>(=60Mpa=8,702psi) |
|                | 5µm: 300bar<br>(=30Mpa=4,351psi)   | 5µm: 300bar<br>(=30Mpa=4,351psi)   | 5µm: 300bar<br>(=30Mpa=4,351psi)   | 5µm: 300bar<br>(=30Mpa=4,351psi)   | 5µm: 300bar<br>(=30Mpa=4,351psi)   |

# Reasons for wide pore columns



*The peak becomes clearer as the pore size increases. Especially at 1000 Å, you can get a clearer peak than the standard 300 Å. And this 1000Å silica gel substrate is rare in the industry.*

Conditions;

Column:

FlexFire C18, 2.6µm (2.0x50mm)

FlexFire WP C4, 2.6µm (2.0x50mm)

FlexFire mAb-RP, 2.6µm (2.0x50mm)

Mobile phase:

A) Water + 0.1%TFA B) Acetonitrile + 0.1%TFA

Gradient:

| min  | mL/min | %A | %B | Curve |
|------|--------|----|----|-------|
| 0.00 | 0.3    | 80 | 20 |       |
| 8.40 | 0.3    | 40 | 60 | 6     |
| 8.42 | 0.3    | 80 | 20 | 6     |

Temperature:

40°C

Detection:

UV280nm

Sample:

1. α-Lactalbumin (0.34mg/mL)

2. β-Lactoglobulin (1.00mg/mL)


Injection volume: 2.0uL




System:

Waters ACQUITY UPLC H-Class PLUS

Mixer:

100uL





### Conditions:

Column: FlexFire WP C4, 2.6µm (2.0x50mm)  
FlexFire AQ C18, 2.6µm (2.0x50mm)

Mobile phase: A) Water + 0.1%TFA  
B) Acetonitrile + 0.1%TFA

### Gradient:

| Time | mL/min | %A | %B | Curve |
|------|--------|----|----|-------|
| 0.00 | 0.3    | 80 | 20 |       |
| 8.40 | 0.3    | 40 | 60 | 6     |
| 8.42 | 0.3    | 80 | 20 | 6     |

Temperature: 70°C

Detection: UV280nm

Sample: 1. Intact Mouse IgG1 (5.0mg/mL)

Injection volume: 1.0µL


System: Waters ACQUITY UPLC H-Class PLUS

Mixer: 100µL

*Wide pore columns are very important for the analysis of compounds with very large molecular weights such as antibodies.*

*Especially for large molecular weight (100KDa-), 1000Å is effective.*

# Separation of catalase (220KDa)



## Conditions:

Column: FlexFire mAb-RP, 2.6µm (2.0x50mm)  
Mobile phase: A) Water + 0.1%TFA B) Acetonitrile + 0.1%TFA

## Gradient:

| min  | mL/min | %A | %B | Curve |
|------|--------|----|----|-------|
| 0.00 | 0.3    | 80 | 20 |       |
| 8.40 | 0.3    | 40 | 60 | 6     |
| 8.42 | 0.3    | 80 | 20 | 6     |

## Temperature:

40°C

## Detection:

UV210nm

## Sample:

1. Ribonuclease A (13.7KDa)
2. Cytochrome C (12.4KDa)
3. Lysozyme (14.3KDa)
4. BSA (66.3KDa)
5. Myoglobin (11.2KDa)
6. Catalase (220KDa)

Injection volume: 2.0µL

## System:

Waters ACQUITY UPLC H-Class PLUS


## Mixer:

100µL

*When catalase with a larger molecular weight is separated and compared with other proteins, the larger the pore size, the better the separation.*

# Separation example by wide pore column

## pBr322 MspI Digest



### Conditions:

Column: FlexFire C18, 1.6 $\mu$ m (2.0x50mm)  
FlexFire WP C18, 1.6 $\mu$ m (2.0x50mm)


Mobile phase: A) 0.1M TEAA, pH7.0 B) 0.1M TEAA, pH7.0/ACN=80/20

### Gradient:

| min  | mL/min | %A   | %B   | Curve |
|------|--------|------|------|-------|
| 0.00 | 0.2    | 42.5 | 57.5 |       |
| 20.0 | 0.2    | 15.5 | 84.5 | 6     |
| 20.1 | 0.2    | 42.5 | 57.5 | 6     |

Temperature: 50°C  
Sample: pBr322 MspI Digest  
Injection volume: 10 $\mu$ L

System: Waters ACQUITY UPLC H-Class PLUS  
Mixer: 100 $\mu$ L



3% Agarose 21

EtBr Staining



*Nucleic acids are not as large in molecular weight as proteins. Therefore, it is possible to analyze even with a normal pore column. However, we have found that using a wide pore column for the analysis of nucleic acids derived from the human body gives very good sensitivity.*

# **Finally**

*Do you notice?*

*Most applications are conditions that can be implemented in LC / MS. We can derive the best analytical conditions for many users.*

*At present, there is a shift from low molecular weight compounds to high molecular weight compounds. We have already developed columns applicable to these compounds. Check the details on the homepage.*

*We cover all the processes ourselves. This is rare worldwide. That's why we can spend time on a lot of research.*

**Get started with FlexFire! !**

