

## 糖類の分析

糖類の分析にはNH<sub>2</sub>カラムを用いることが多く、分離はHILICと同じ挙動を示します。そこで水の多い移動相においても耐久性のある「Develosil ANIDIUS」を用いた糖類の分析例を紹介します。




Fig.1 Analysis of Saccharides

### Conditions;

Column: Develosil ANIDIUS  
(4.6x250mm)

Mobile phase: Acetonitrile/Water=80/20

Flow rate: 1.0ml/min

Temperature: 40°C

Detection: RI

Sample: Fructose, Glucose, Sucrose,  
Maltose

Injection volume: 20uL

System: 日本分光製




Fig.2 Analysis of Maltooligosaccharides

### Conditions;

Column: Develosil ANIDIUS  
(4.6x250mm)

Mobile phase: Acetonitrile/Water=65/35

Flow rate: 1.0ml/min

Temperature: 40°C

Detection: RI

Sample: Glucose, Maltose, Maltotriose,  
Maltotetraose, Maltopentaose  
maltohexaose, Maltoheptaose

Injection volume: 20uL

System: 日本分光製