



No.097



Column: Develosil HSR AQ C18, 5um (4.6x150mm)
Develosil HSR AQ C18, 2um (2.0x50mm)
Mobile phase: A) Water + 0.1%HCOOH
B) Acetonitrile + 0.1%HCOOH
Flow rate: 1.0mL/min (5um)
0.474mL/min (2um)
Temperature: 40°C
Detection: 275nm
Sample: 1.Sulfadiazine (0.08mg/mL)
2.Sulfathiazole (0.10mg/mL)
3.Sulfapyridine (0.10mg/mL)
4.Sulfamerazine (0.09mg/mL)
5.Sulfamethazine (0.09mg/mL)
6.Sulfamethoxy pyridazine (0.09mg/mL)
7.Sulfachloropyridazine (0.10mg/mL)