


## 抗がん剤(1)



### Conditions:

Column: Develosil XG-C1-5, XG-C18LC-5, XG-C18M-5, XG-C30M-5

Column size: 4.6x150mm

Mobile Phase: 0.2%Phosphoric acid/Methanol=95/5

Flow rate: 1.0mL/min

Temperature: 40°C

Sample: 1:Potassium Oxonate (0.17mg/mL)

2:5-Fluorouracil (0.17mg/mL)

3:Gimeracil (0.17mg/mL)

4:5-Fluoro-1-(tetrahydro-2-furfuryl)uracil (0.15mg/mL)

Injection Volume: 2.0μL

Detection: PDA (Max Abs.)