


# Develosil Application Data

## No.066 Analgestic drug/鎮痛薬



### Analytical Conditions:

Column: Develosil HSR C18,5um (4.6x150mm)

Mobile phase: Acetonitrile/25mM NH<sub>4</sub>H<sub>2</sub>PO<sub>4</sub>, pH3.5=15/85

Flow rate: 1.0mL/min

Temperature: 40°C

Detection: UV240nm

Sample: 1.Uracil (5.0ug/mL)

2.Acetaminophene (4.2ug/mL)

3.Caffeine (0.64mg/mL)

4.Aspirine (0.59mg/mL)

Injection volume: 2.0uL

System: Waters alliance